





CABLE 2m STANDARD

4 x Ø 3.30 ± 15  
M4 - 6H ± 12  
✓ Ø 4.20 X 90°, Near Side  
MAX SCREW TORQUE 0.6N.m

4 x Ø 5.40 THRU



AXIS OF MOTION  
(IF MOUNTED VERTICALLY  
THIS SHOULD BE UPWARDS)

|                                                        |                                                     |              |      |   |        |          |
|--------------------------------------------------------|-----------------------------------------------------|--------------|------|---|--------|----------|
| TITLE                                                  | INSTALLATION OF NPS-Z-250B / NPS-X-250B PIEZO STAGE |              |      |   | APRVD. |          |
| DRG. No.                                               | QG-0534-IN-100002                                   | SHEET 1 OF 1 | ISS. | A | DRAWN  | J.C.     |
| MATL.                                                  | ALUMINIUM STAGE BODY                                |              |      |   | SCALE  | 1:2      |
| FINISH                                                 | NICKEL                                              |              |      |   | DATE   | 14.05.25 |
| PRIOR Scientific Instruments Ltd. Fulbourn, Cambridge. |                                                     |              |      |   |        |          |

|                                             |  |                                 |
|---------------------------------------------|--|---------------------------------|
| 3RD ANGLE                                   |  | DO NOT SCALE<br>IF IN DOUBT ASK |
| TOLERANCE UNLESS OTHERWISE STATED<br>±0.1mm |  |                                 |
| BREAK ALL SHARP EDGES                       |  |                                 |

| ECN     | REV | CREATED | APP | DATE | DETAILS     |
|---------|-----|---------|-----|------|-------------|
| 1st Rel | A   | J.C.    |     |      | 1st Release |
|         |     |         |     |      |             |
|         |     |         |     |      |             |
|         |     |         |     |      |             |
|         |     |         |     |      |             |