



NOTES

WHEN USING THE STAGE LIFTING AGAINST GRAVITY THE STAGE MUST BE MOUNTED WITH THE ARROW UPWARDS TO MAINTAIN PRESSURE ON THE PIEZO ACTUATOR WITHIN.

TITLE	XP200 200um	APRVD.		3RD ANGLE 	DO NOT SCALE IF IN DOUBT AS
DRG. No.	QGXP200	SHEET 1 OF 1	ISS. A	DRAWN J.C.	
MATL.			SCALE	2:1	TOLERANCE UNLESS OTHERWISE STATED $\pm 0.1\text{mm}$
FINISH			DATE	03.02.25	
PRIOR Scientific Instruments Ltd. Fulbourn, Cambridge.					BREAK ALL SHARP EDGES

ECN	REV	CREATED	APP	DATE	DETAILS
1st Rel	A	J.C.			1st Release